

LabSystems

LABORATORY SYSTEMS GROUP PTY LTD


THE ACADEMY SCHOOL FUME CUPBOARD

*LabSystems Group Introduces "The Academy Education Series"
Bringing your Secondary Schools or Colleges,
Science and Language Centers Safety with Style into the 21st Century.*

The Academy Series ducted fume cupboard- This Academy Education Series School type fume cupboard has been redesigned specifically for the secondary school and college market. It is available in both single and rear loading (double sided) configurations and uses the most modern polymer and material technology in its construction and complies fully with Australian Standards AS 2243.8-2006.

The Academy Features;

- Modern Ergonomic Design
- Robust construction
- Chemical resistant work surface and main chamber
- Angled glass sash for ergonomic viewing.
- Active bypass
- Academy Safety Control System with pre and post purge as per AS/NZ 2243.8 - 2009
- Available in 915mm and 1200mm width
- Single sided or Double sided (rear loading) Prep room model.



Please don't hesitate to contact our sales office on (03) 8720 9000, if you require any further information or view our website for our complete range of products and services.



WORLD RECOGNISED
ACCREDITATION

JAS-ANZ



Laboratory Systems Group Pty Ltd

ABN 35 796 401 519

Unit 6, 144-150 Canterbury Road, Kilsyth, VIC 3137

Telephone: 03 8720 9000 Fax: 9761 7350

Email: sales@labsystemsgroup.com.au website: www.labsystemgroup.com.au