

# POLYSTORE

## CHEMICAL STORAGE CABINETS

"The Polystore chemical storage cabinets offer a range of cabinet sizes for all corrosive storage requirements"

### Features:

- Fire retardant Polypropylene construction (tested to AS 1530.3-1989)
- Liquid tight sump
- No internal metal parts means no corrosion that is associated with metal cabinets
- Adjustable and removable shelves
- 50mm diameter vents
- Self-closing doors
- Lockable doors
- Comply with AS 3780

### Model Specification:

| CAT NO.     | Volume (litre) | Height (mm) | Width (mm) | Depth (mm) |
|-------------|----------------|-------------|------------|------------|
| CSCO50L-UB  | 50             | 750         | 600        | 450        |
| CSCO100L-UB | 100            | 800         | 600        | 450        |
| CSCO150L-UB | 150            | 1000        | 1025       | 500        |
| CSCO200L-UB | 200            | 1500        | 1025       | 550        |

Note: FS =Free Standing  
UB = Underbench - D = Drum Storage



Note: Specific size requirements can be made to order.

