

LABSTORE

CHEMICAL STORAGE CABINETS

The LabStore cabinets offer a range of cabinet sizes for hazardous liquid storage requirements.

Features:

- Manufactured from 0.9mm powder-coated zincanneal
- Liquid tight sump
- Self-closing doors
- Adjustable shelves
- Lockable doors
- 50mm diameter vents
- Double walled with 40mm insulating air space
- All cabinets comply with AS 1940, AS 3780, AS 4326, AS 4452 and AS 2714
- Lockable doors
- All cabinets comply with AS 3780

Model Specification:

CAT NO.	Volume (litre)	Height (mm)	Width (mm)	Depth (mm)
AU25714	30	800	595	465
AU52712	60	1005	595	465
AU25302	160	1230	1100	465
AU25452	250	1760	1100	465
AU25602	250	1760	870	870

AU25452 Volume - 9 x 20 Litre Drums
AU25602 Volume - 12 x 20 Litre Drums

Suitable for Storage of:

- Flammable Substances
- Corrosive Substances
- Organic Peroxides
- Oxidizing agents
- Toxic Compounds



Note: Specific size requirements can be made to order.

