

LABORATORY SYSTEMS GROUP DOWNDRAFT DISSECTION WORKBENCH

Model # LSG-DFB 1500PP

Rigid Industrial Polypropylene Construction

Dimensions: 1500mm wide x 940mm high x 750mm deep (Splash back to 1800mm)

Total downdraft worktop, inner chamber has drain plenum

Inner chamber totally accessible for cleaning

Large polypropylene wash sink 300mm x 300mm x 200mm deep (LH or RH)

Removable downdraft worktops

Front access for duct connection and plumbing

Standard splash back

Side or rear exhaust duct entry

OPTIONS

Hot and Cold water fittings with single lever mixing valve

Hand held spray wash down fitting

Custom sizes: Custom size can be manufactured on request

