

AirClean® Systems

"THE FUME CONTROL EXPERTS"

POLYPROPYLENE - CLASS 100 HORIZONTAL CLEAN BENCH

AirClean® Systems Horizontal Laminar Flow Clean Benches are constructed from all-white, structural polypropylene. This provides a seamless design and exceptional Class 100 laminar flow air in a horizontal direction. Uniform air exists from the rear of the hood where it is filtered and then flows forward to the face of the hood. Built-in fluorescent lighting is provided to illuminate the work surface. Units are available in standard 1320mm, 1524mm, 1829mm and 2590mm widths while custom widths and depths are available.



Features:

- All polypropylene construction - NO MORE RUST!
- Seamless, smooth construction minimizes air turbulence
- HEPA clean air
- Minihelic gauge continuously monitors filter condition
- Optional cup sink
- Optional variable air flow control
- Built-in fluorescent lighting
- Available in 110V or 220V
- Custom sturdy carts
- Meets FS209E Class 100
- All-white surfaces are easy to clean

MODEL	INSIDE DIMENSIONS	OUTSIDE DIMENSIONS
AC4000HLF	119mm x 502mm x 635mm	1320mm x 762mm x 1016mm
AC5000HLF	1473mm x 559mm x 635mm	1524mm x 762mm x 1016mm
AC6000HLF	1778mm x 559mm x 635mm	1829mm x 762mm x 1016mm
AC8000HLF	2464mm x 584mm x 635mm	2590mm x 762mm x 1016mm

