

## IsoGuard

### Isolation Containment Chamber





The **IsoGuard** Isolation Chamber provides the perfect processing micro-environment for a host of containment and isolation applications. This unit is fully self contained and suitable for PC-2,3,4, QC, 2,3,4 Cleanrooms, toxic powders and Bio containment.

#### Features

- GUARD™ Electronic control system and HEPA monitoring
- Lockable pass-thru chamber
- Stainless steel or structural Polypropylene construction
- HEPA or ULPA filtration
- Optional carbon filtration inlet and outlet
- Bag in/bag out HEPA changing
- Alternative glove materials available
- Fluorescent light
- Complete front opening access for initial loading
- Polypropylene model has exceptional resistance to aggressive decontamination media

# IsoGuard

## Isolation Containment Chamber



### Optional Features

- Trolley
- UV Light for decontamination
- Decontamination Kits
- Custom sizes, access & fittings
- Complete metal free interior
- Thimble style exhaust connection (total exhaust)
- Gases & GPO's

