


# ChemGuard™ Filtered Filtered Storage Cabinet





Concerned about chemical smells in your lab? The Laboratory Systems Group ChemGuard™ filtered storage cabinet is the perfect solution—and requires no ducting. Our vented, filtered storage cabinet is available in a standard configuration to suit most chemical storage needs, and offers adjustable shelving and high chemical resistance polypropylene design. (custom sizes on application )The ChemGuard™ range have fully filtered and monitored, recirculating exhaust with no ductwork required, and a range of activated carbon filters are available to suit most chemicals. To discuss your chemical storage needs and find out whether the ChemGuard™ is right for you, contact us today.

## Features

- Fully filtered & recirculating exhaust
- No ducting required
- Chemical resistant polypropylene design
- Adjustable Bunded PP shelving
- Range of activated carbon filters available
- Clear viewing doors for easy access & selection of goods
- Chemical monitoring of filter saturation
- Not suitable for all chemical classes—contact us to find out more
- LED strip lighting
- Door lock

Rev 17/06

# ChemGuard™ Filtered Filtered Storage Cabinet



MODEL	VOLUME	HEIGHT	WIDTH	DEPTH
CG9521RC		2100MM	900MM	500MM

Laboratory Systems Group Pty Ltd  
24 Melba Avenue, Lilydale, Victoria, 3140 Australia

Phone: +61 1300 LAB SYS  
1300 522 797  
Fax: +61 3 9738 7741  
Email: sales@labsystemsgroup.com.au  
Web: www.labsystemsgroup.com.au

