

ChemGuard™ Static

Static Storage Cabinet



The Laboratory Systems Group ChemGuard™ Static Storage Cabinet provides a static sealed enclosure for hazardous samples. With a sleek, appealing design, the ChemGuard™ static storage cabinet is manufactured using strong, fire retardant polypropylene and features toughened glass doors. Should you require an alternate material or size, our storage cabinet is fully customisable and can be tailor-made to suit your requirements. The ChemGuard™ static is rigorously tested to Australian Standard AS/NZS 1530.3 for material flammability, AS/NSW 2243.10 for sump capacity and AS/NZS 1668.2 for total exhaust models.


Features

- Aesthetically appealing design
- Static sealed enclosure
- Fire retardant polypropylene construction
- Liquid tight sump
- Adjustable & removable shelves
- Toughened glass doors
- Custom sizes & materials available

Rev 17/06

ChemGuard™ Static

Static Storage Cabinet



| MODEL | HEIGHT | WIDTH | DEPTH |
|---------|--------|-------|-------|
| CG9521S | 2100MM | 900MM | 500MM |