

PowderGuard™

Ductless Fume Cupboard



The PowderGuard™ ductless fume cupboard is the perfect solution for the handling of low-toxicity powders, substances and chemicals within a safe, contained environment. Ideal for weighing and manipulation of pharmaceutical powders, the PowderGuard™ protects the operator from hazardous compounds using an advanced HEPA filtration system. The PowderGuard™ requires no ducting, allowing for portability, and features a state-of-the-art control system that alerts user to low airflow and filter replacement needs.

The PowderGuard™ is designed and manufactured in Australia, and is available in a variety of sizes to suit your needs. A number of optional features are available, including the addition of carbon filtration. The PowderGuard™ recirculating cabinet has been rigorously tested to ensure the safety of you and your laboratory.


Features

- Filtered, ductless design
- Designed & made in Australia
- Compact & modern for smaller environments
- Minimal installation
- Low airflow & filter replacement alarm
- Features HEPA filtration system
- Safe & simple pre-filter replacement process
- LED lighting

Rev 17/06

PowderGuard™

Ductless Fume Cupboard



Optional Features

- Polypropylene sink & drain
- Gas, water & GPO services (with base option—as shown)
- Custom base inserts
- Mobile under bench
- Electric height-adjustable bench (shown)

DIMENSIONS	PowderGuard 800™	PowderGuard 1200™
HEIGHT—EXTERNAL	1000mm	1000mm
WIDTH—EXTERNAL	805mm	1205mm
DEPTH—EXTERNAL	600mm	600mm