

ChemGuard™ Vented Vented Storage Cabinet



Our ChemGuard™ Vented Storage Cabinet provides a perfect ventilated storage solution for odourous samples. Available in a wide range of cabinet sizes, the Vented Storage Cabinet allows for safe sample storage and is ideal for chemistry, pathology, pharmaceutical and biological laboratories.

With a sleek, appealing design, the Vented Storage Cabinet is manufactured using strong, fire retardant polypropylene and features toughened glass doors. Should you require an alternate material or size, our storage cabinet is fully customisable and can be tailor-made to suit your requirements. The Vented Storage Cabinet is rigorously tested to Australian Standard AS/NZS 1530.3 for material flammability, AS/NSW 2243.10 for sump capacity and AS/NZS 1668.2 for total exhaust models.


Features

- Aesthetically appealing design
- Fire retardant polypropylene construction
- Liquid tight sump
- Adjustable & removable shelves
- Exhaust spigot to suit many sizes
- Toughened glass doors
- Custom sizes & materials available
- Available in filtered/recirculating variation with bonded carbon filters for chemical & odour removal

Rev 19/08

ChemGuard™ Vented

Vented Storage Cabinet



MODEL	HEIGHT	WIDTH	DEPTH
CG9521V	2100MM	900MM	500MM