

Drying Racks

Glassware Draining Racks



Laboratory Systems Group draining racks are made from a sturdy (PVC) coated steel frame, with ample room for hanging and drying a range of glassware and test tube sizes. Available in several standard sizes or able to be custom made for your laboratory. The clear acrylic back is easily fixed to the rear of your sink, or any substrate. Featuring a 45° lip at the base, water is directed away from your wall.



Features

- Available in several standard sizes :
 - 600 – one rack
 - 1200 – two racks
 - 1500 – two racks
 - 1800 – three racks
- PVC coated steel glassware support spigots
- Acrylic back with diversion lip at base
- Pre-drilled fixing holes
- Sturdy plastic coated steel frame

Rev 18/04

Drying Racks

Glassware Draining



Options:

- Back board can be manufactured from materials other than acrylic
- Stainless steel frame
- Custom shaped backboard to avoid existing tapware
- Custom widths to suit your available space

MODEL	WIDTH (mm)	HEIGHT (mm)	FRAMES	WEIGHT (kg)
ADR-6	600	930	1	3
ADR-12	1200	930	2	6
ADR-15	1500	930	2	8
ADR-18	1800	930	3	10