

Polystore

Corrosive Storage Cabinet



The Polystore corrosive storage cabinet is the perfect solution for containing corrosive chemicals and substances in your laboratory. Available in a wide range of cabinet sizes, the Polystore offers a safe, stylish storage option featuring no internal metal parts, fire retardant polypropylene construction, self closing doors and liquid tight sump. The Polystore corrosive storage cabinet has undergone many years of development and testing, and is fully compliant to Australian Standard AS/NZS 3780 and AS/NZS 2243.10.

Our Polystore corrosive storage cabinets can be custom made to your requirements and dimensions, and are available with adjustable and removable shelving, and self-closing and lockable doors. If you are looking for a chemical storage cabinet to keep you, your laboratory and your staff safe, contact one of our Team to find out more about Polystore today.


Features

- Multiple sizes available
- Fire retardant polypropylene construction
- Liquid tight sump
- No internal metal parts
- Stylish, ergonomic design
- Available with right or left-hand opening
- Self-closing & lockable door
- Fully compliant to AS/NZS 3780 & AS/NZS 2243.10



Polystore

Corrosive Storage Cabinet



Optional Features

- Adjustable and removable shelves
- Combination acid/base units in the one cabinet using divider—available 100, 150, 200 & 250 litre models
- Fully customised to laboratory requirements
- Mechanical ventilation knock-out

MODEL	VOLUME	HEIGHT	WIDTH	DEPTH	WEIGHT	PART NUMBER	EXTRA SHELF PART #
CSC 50 LITRE Hinged on RHS	50 LITRES	720mm	600mm	400mm	19kg	101489	101514
CSC 50 LITRE Hinged on LHS	50 LITRES	720mm	600mm	400mm	19kg	101578	101514
CSC 100 LITRE	100 LITRES	720mm	935mm	400mm	29kg	101613	101621
CSC 100 LITRE With Divider	100 LITRES	720mm	935mm	400mm	33kg	110535	124798
CSC 150 LITRE	150 LITRES	1000mm	1000mm	500mm		115705	115701
CSC 150 LITRE With Divider	150 LITRES	1000mm	1000mm	500mm		118306	118313
CSC 200 LITRE	200 LITRES	1250mm	1000mm	450mm	87kg	102285	102293
CSC 200 LITRE With Divider	200 LITRES	1250mm	1000mm	450mm		108560	108563
CSC 250 LITRE	250 LITRES	1760mm	1100mm	450mm		108341	102293
CSC 250 LITRE With Divider	250 LITRES	1760mm	1100mm	450mm		116182	